



# Anti-DPP4 (C-terminal) polyclonal antibody (DPAB-DC4189)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Immunogen	A synthetic peptide corresponding to C-terminus of zebrafish DPP4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Zebrafish
Conjugate	Unconjugated
Applications	WB, ELISA,
Format	Liquid
Size	100 µg
Buffer	In PBS, pH 7.2 (50% glycerol, 0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

Keywords	DPP4
----------	------

## GENE INFORMATION

Entrez Gene ID [1803](#)

---

UniProt ID [P27487](#)

---