



Anti-FABP polyclonal antibody (DPAB-DC4228)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Target	FABP
Immunogen	Native purified fatty acid binding protein from human heart tissue.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB, IP, EIA,
Procedure	Cardiac markers Antibodies
Format	Liquid
Size	150 µg
Buffer	In PBS, pH 7.5 (50% glycerol, 0.01% thimerosal)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	FABP
-----------------	------