



Anti-Fibrinogen polyclonal antibody (DPAB-DC4239)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	This product is suitable for use as a precipitating antibody. Cross-reactivity to other species may exist. Antibody activity to human fibronectin or serum proteins is not present.
Target	Fibrinogen
Immunogen	Human FG.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB, ELISA, ID,
Procedure	Cardiac markers Antibodies
Format	Liquid
Size	2 ml
Buffer	In 20 mM sodium phosphate buffer, 140 mM NaCl, pH 7.3 (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords

FG

GENE INFORMATION

Entrez Gene ID

[2243](#)

UniProt ID

[P02671](#)
