



# Anti-MSH-alpha polyclonal antibody (DPAB-DC4548)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Specificity</b>	This antibody may also react with Mouse, Rat and Rabbit
<b>Immunogen</b>	A synthetic peptide (conjugated with BSA) corresponding to MSH-alpha.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	IHC-P, IHC-Fr,
<b>Format</b>	Liquid
<b>Size</b>	50 µl
<b>Buffer</b>	In Tris-HCl buffer (1% BSA, 0.09% sodium azide)
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

## BACKGROUND

<b>Introduction</b>	α-Melanocyte-stimulating hormone (alpha-MSH; α-MSH), also known as alpha-melanotropin, alpha-melanocortin, or alpha-intermedin, is a naturally occurring endogenous peptide hormone of the melanocortin family, with a tridecapeptide structure and the amino a
---------------------	---

**Keywords**

$\alpha$ -MSH; MSH-alpha

---