



Anti-PSD-95 (N-terminal) polyclonal antibody (DPAB-DC4667)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to N-terminus rat PSD-95.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Bovine, Human, Mouse, Primates, Rat, Zebrafish
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB (Tissue lysate),
Format	Liquid
Size	100 µl
Buffer	In 10 mM HEPES, 150 mM NaCl, pH 7.5 (50% glycerol, 10% BSA)
Preservative	None
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Keywords	PSD-95
-----------------	--------

GENE INFORMATION

Entrez Gene ID [1742](#)

UniProt ID [G5E939](#)
