



Anti-Thyroid Stimulating Hormone polyclonal antibody (DPAB-DC4824)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	human Alpha TSH, Beta TSH, and Intact TSH.
Target	Thyroid Stimulating Hormone
Immunogen	Purified human TSH.
Source/Host	Goat
Species Reactivity	Human
Conjugate	Unconjugated
Applications	RIA,
Format	Liquid
Size	1 ml
Buffer	In buffer containing 0.1% sodium azide
Preservative	0.1% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Thyroid-stimulating hormone (also known as thyrotropin, TSH, or hTSH for human TSH) is a hormone that stimulates the thyroid gland to produce thyroxine (T4), and then triiodothyronine (T3) which stimulates the metabolism of almost every tissue in the body
---------------------	---

KeywordsTSH; Thyroid-stimulating hormone

GENE INFORMATION**Entrez Gene ID**[12640](#)

UniProt ID[A2AVN3](#)
