



Anti-LRIG3 (aa 1020-1119) polyclonal antibody (DPAB-DC490)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	LRIG3 (leucine-rich repeats and immunoglobulin-like domains 3) is a protein-coding gene. Diseases associated with LRIG3 include ependymoma, and synovitis. An important paralog of this gene is NYX.
Immunogen	LRIG3 (NP_700356, 1020 a.a. ~ 1119 a.a) partial recombinant protein with GST tag. The sequence is LDFSANPEPASVASSNSFMGTFGKALRRPHLDAYSSFGQPSCQPRAFYLKAHSSPDLDS GSEEDGKERTDFQEENHICTFKQTLENYRTPNPFQSYDLDT
Source/Host	Mouse
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB (Recombinant protein), ELISA,
Size	50 µl
Buffer	50 % glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [LRIG3 leucine-rich repeats and immunoglobulin-like domains 3 \[Homo sapiens \(human\) \]](#)

Official Symbol	LRIG3
Synonyms	LRIG3; leucine-rich repeats and immunoglobulin-like domains 3; LIG3; leucine-rich repeats and immunoglobulin-like domains protein 3; LIG-3;
Entrez Gene ID	121227
Protein Refseq	NP_001129523
UniProt ID	Q6UXM1
Chromosome Location	12q14.1