



Anti-GAREML (C-terminal) polyclonal antibody (DPAB-L20436)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	FAM59B antibody was raised against a 19 amino acid peptide from near the carboxy terminus of human FAM59B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	immunohistochemistry: suitable
Format	Buffered aqueous solution
Size	100 µg
Buffer	Supplied in PBS with 0.02% sodium azide.
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C

GENE INFORMATION

Gene Name	GAREML GRB2 associated, regulator of MAPK1-like [Homo sapiens (human)]
Official Symbol	GAREML
Synonyms	GAREML; GRB2 associated, regulator of MAPK1-like; FAM59B; GRB2-associated and

regulator of MAPK protein-like; protein FAM59B; GRB2-associated and regulator of MAPK1-like; family with sequence similarity 59, member B

Entrez Gene ID	150946
-----------------------	------------------------

Protein Refseq	NP_001161713
-----------------------	------------------------------

UniProt ID	Q75VX8
-------------------	------------------------
