



Anti-Dengue Virus Type 1 polyclonal antibody (DPAB0110)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit polyclonal to Dengue Virus 1+2+3+4
Specificity	Reacts with Dengue Virus Types 1, 2, 3 and 4.
Target	Dengue Virus Type 1,2,3,4
Immunogen	Viral particles from C6/36 [mosquito cell line] amplified dengue type 2 virus.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Ammonium sulphate precipitated.
Conjugate	Unconjugated
Applications	ICC/IF
Positive Control	dengue type III infected BHK21 cells.
Format	Liquid
Concentration	7.34mg/ml
Size	250 µl
Buffer	PBS, pH 7.2
Preservative	0.1% Sodium Azide

Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.
----------------	--

BACKGROUND

Introduction	Dengue and dengue hemorrhagic fever (DHF) are caused by one of four closely related, but antigenically distinct, virus serotypes (DEN1, DEN2, DEN3, and DEN4), of the genus Flavivirus.
Keywords	Dengue Virus 1+2+3+4;Dengue Virus Type 1,2,3,4;Group IV;Flaviviridae;Flavivirus;Dengue virus;Rabbit Anti-Dengue Virus Type 1,2,3,4 Polyclonal Antibody;Rabbit Anti-Dengue Virus Type 1,2,3,4;Polyclonal Antibody;Rabbit Anti-Dengue Virus Type 1,2,3,4 Polyclonal;Antibody;Rabbit Anti-Dengue Virus Type 1,2,3,4 Antibody;Polyclonal;Rabbit Anti-Dengue Virus 1+2+3+4 Antibody;Rabbit Anti-Dengue Virus 1+2+3+4 Polyclonal;Anti-Dengue Virus 1+2+3+4 Antibody;Anti-Dengue Virus 1+2+3+4 Polyclonal;
