



Rabbit Anti-HBsAg PreS1 (subtypes ad & ay) Polyclonal antibody (DPAB0225)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Monospecific, reacts only with Hepatitis B surface antigen including the pre-S1 epitope. Non-reactive with normal human serum.
Target	HBV Surface Antigen
Immunogen	Hepatitis B surface antigen purified from human serum. Mixture of subtypes ad & ay.
Source/Host	Rabbit
Species Reactivity	HBV
Purification	IgG fraction covalently coupled with the N-Hydroxysuccinimide ester of biotin under mild conditions to give a high degree of substitution
Conjugate	Biotin
Applications	Suitable for use with avidin and streptavidin amplification systems for ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Concentration	4-5mg/ml (OD280nm, E0.1% = 1.4)
Size	1 ml
Buffer	0.01M PBS, pH 7.2 This product contains no stabilizing proteins.
Preservative	0.1% Sodium Azide
Storage	Short-term (up to 6 months) store at 2–8°C. Long term, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Introduction

Hepatitis B Virus (HBV) infection induces a disease state characterised by liver damage, inflammation and viral persistence. Infection also increases the risk of hepatocellular carcinoma. HBV belongs to the Hepadnaviridae family of viruses. Its genome consists of partially double stranded circular DNA. The DNA is enclosed in a nucleocapsid, or core antigen (HBcAg), which is surrounded by a spherical envelope (surface antigen or HBsAg). The core antigen shares its sequences with the e antigen (HBeAg) but no cross reactivity between the two proteins has been observed. The HBV genome also encodes a DNA polymerase that also acts as a reverse transcriptase.

Keywords

HBsAg; HBV major surface antigen; HBV surface antigen; Hepatitis B Virus major surface antigen; Major surface antigen; S; Hepatitis B Surface Antigen; Hepatitis B Virus Surface Antigen; Hepatitis B virus; Orthohepadnavirus; Hepadnaviridae