



Sheep anti-Proc Polyclonal antibody (DPAB21455H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Proc Polyclonal Antibody
Target	Proc
Source/Host	Sheep
Conjugate	Unconjugated
Applications	ELISA
Molecular Weight	150000
Concentration	22.3 mg/mL
Buffer	50% Glycerol/H ₂ O (vol/vol)
Preservative	None
Storage	-20 °C

GENE INFORMATION

Gene Name	Proc protein C [Mus musculus]
Official Symbol	PROC
Synonyms	PC
Entrez Gene ID	19123

Protein Refseq	NP_032960.3
UniProt ID	P33587
Chromosome Location	18; 18 B3
Pathway	Cell surface interactions at the vascular wall, organism-specific biosystem;Common Pathway, organism-specific biosystem;Complement and Coagulation Cascades, organism-specific biosystem;Complement and coagulation cascades, organism-specific biosystem;Compl
Function	calcium ion binding;catalytic activity;hydrolase activity;peptidase activity;serine-type endopeptidase activity;