



Anti-Robenidine polyclonal antibody (DPAB27943)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Robenidine Polyclonal Antibody
Target	Robenidine
Immunogen	Robenidine-BTG
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	IA
Format	Liquid/Ig Fraction
Concentration	13.41mg/ml
Size	1 mg
Buffer	20mM Phosphate, 150mM Sodium Chloride
Preservative	0.09% Sodium Azide
Storage	Store it under sterile conditions at -20°C upon receiving. Recommend to pack the protein into smaller quantities for optimal storage.

BACKGROUND

Introduction	Robenidine is a coccidiostat.
---------------------	-------------------------------

Keywords

1,2-Bis[(4-chlorophenyl)methyleneamino]guanidine ; Sheep Anti-Robenidine Polyclonal Antibody; Anti-Robenidine Polyclonal Antibody; Robenidine Polyclonal Antibody Sheep Anti-Robenidine PAb; Anti-Robenidine PAb; Robenidine PAb; Sheep Anti-Robenidine Antibod
