



Anti-Spectinomycin polyclonal antibody (DPAB27959)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Spectinomycin Polyclonal Antibody
Target	Spectinomycin
Immunogen	Spectinomycin-BSA
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	IA
Format	Liquid/Ig Fraction
Concentration	12.36mg/ml
Size	1 mg
Buffer	20mM Phosphate, 150mM Sodium Chloride
Preservative	0.09% Sodium Azide
Storage	Store it under sterile conditions at -20°C upon receiving. Recommend to pack the protein into smaller quantities for optimal storage.

BACKGROUND

Introduction Spectinomycin is an aminocyclitol antibiotic, closely related to the aminoglycosides, produced

by the bacterium *Streptomyces spectabilis*. Aminoglycosides are a family of bacterial antibiotics that are used in the treatment of specific bacterial infections. They display a concentration dependent killing action and are active against a wide range of aerobic Gram-negative bacilli. Aminoglycosides are molecules that are comprised of an amino group and a sugar group. They operate by inhibiting the bacteria from producing proteins vital to its growth. More specifically, they bind to the bacterial 30S ribosomal subunit where they prevent the translocation of the peptidyl-tRNA from the A-site to the P-site, subsequently giving rise to a misreading of mRNA resulting in the inhibition of protein synthesis. This consequently results in a disruption to the integrity of the bacterial cell membrane. In addition to their use to prevent bacterial infection, aminoglycosides have been used as growth promoters in food producing animals.

Keywords

spectinomycin hydrochloride; Sheep Anti-Spectinomycin Polyclonal Antibody; Anti-Spectinomycin Polyclonal Antibody; Spectinomycin Polyclonal Antibody Sheep Anti-Spectinomycin PAb; Anti-Spectinomycin PAb; Spectinomycin PAb; Sheep Anti-Spectinomycin Antibody
