



Anti-Estriol-6 polyclonal antibody (DPAB3472)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal antibody to Human Estriol-6 antigen.
Immunogen	Estriol-6-BSA
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Format	Neat serum.
Preservative	None
Storage	4-8°C short term; -20°C long term.

BACKGROUND

Introduction	Estriol is an estrogenic hormone, a metabolite of estradiol and usually the predominant estrogenic metabolite found in the urine of pregnant women.
Keywords	oestriol; Estriol; Estriol-6