



Anti-Sulfaguanidine polyclonal antibody (DPAB4213)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit polyclonal antibody to sulfaguanidine.
Immunogen	Sulfaguanidine-KLH Conjugate
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Buffer	0.01 mol/L pH 7.4 Phosphate Buffer Solution.
Preservative	None
Storage	Short-term store at 2-8°C, Long-term store at -20°C.

BACKGROUND

Introduction	Sulfaguanidine is a sulfonamide.
Keywords	p-Aminobenzenesulfoguanidide; Resulfon; RP 2275; rp2275; Ruocid; SG; S-Guanidine; Shigatox; Suganyl; Sulfaguamide; Sulfaguanidin; Sulfaguanil; Sulfaguine; Sulfanilamide, N(sup1)-amidino-; Sulfanilamide, N1-(diaminomethylene)-; Sulfanilamide, N1-amidino-; Sulfanilguanidine; sulfanilyl-guanidin; Sulfentidine; Sulfoguanidine; Sulfoguanil; Sulfoguanyl; Sulfoguenil; Sulfoquanidine; Sulgin