



Anti-ARPC5L (aa 1-136) polyclonal antibody (DPABH-00001)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein fragment corresponding to a region within amino acids 1 and 136 of Human ARPC5L (UniProt Q9BPX5).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat, Human
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	50 µl
Buffer	pH: 7,00; Constituents: 0.75% Glycine, 1.21% Tris, 20% Glycerol
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	May function as component of the Arp2/3 complex which is involved in regulation of actin polymerization and together with an activating nucleation-promoting factor (NPF) mediates the formation of branched actin networks.
--------------	--

KeywordsARPC5L; actin related protein 2/3 complex, subunit 5-like

GENE INFORMATION**Entrez Gene ID**[81873](#)

UniProt ID[Q9BPX5](#)
