



Anti-CCDC61 (aa 43-73) polyclonal antibody (DPABH-00442)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The specific function of CCDC61 remains unknown.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within N terminal amino acids 43-73 of Human CCDC61
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	1 ea
Buffer	Constituent: PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Gene Name [CCDC61 coiled-coil domain containing 62 \[Homo sapiens \]](#)

Official Symbol	CCDC61
Synonyms	CCDC61; coiled-coil domain containing 61; coiled-coil domain-containing protein 61;
Entrez Gene ID	729440
Protein Refseq	NP_001254652.1
UniProt ID	Q9Y6R9