



Anti-SLC35F5 (aa 186-195) polyclonal antibody (DPABH-00827)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Putative solute transporter.
Immunogen	Synthetic peptide: C-EKTPKKSRVR, corresponding to internal sequence amino acids 186-195 of Human SLC35F5 (NP_079457.2)
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	200 µl
Buffer	pH: 7.30; Constituents: Tris buffered saline, 0.5% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [SLC35F5 solute carrier family 35, member F6 \[Homo sapiens \]](#)

Official Symbol	SLC35F5
Synonyms	SLC35F5; solute carrier family 35, member F5; solute carrier family 35 member F5; HCV NS5A-transactivated protein 3; hepatitis C virus NS5A-transactivated protein 3;
Entrez Gene ID	80255
Protein Refseq	NP_079457.2
UniProt ID	A0A024RAD4