



Anti-OR3A2 (C-terminal) polyclonal antibody (DPABH-00927)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Odorant receptor.
Immunogen	Synthetic peptide derived from C terminal of Human OR3A2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituents: 49% PBS, 50% Glycerol, 0.88% Sodium chloride. Note: PBS without Mg2+ and Ca2+
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C.

GENE INFORMATION

Gene Name	OR3A2 olfactory receptor, family 3, subfamily A, member 3 [Homo sapiens]
-----------	--------------------------------------------------------------------------------------------

Official Symbol	OR3A2
Synonyms	OR3A2; olfactory receptor, family 3, subfamily A, member 2; OR228; OLFRA04; OR17-14; OR17-228; olfactory receptor 3A2; olfactory receptor 17-228; olfactory receptor OR17-14;
Entrez Gene ID	4995
Protein Refseq	NP_002542.3
UniProt ID	P47893
Pathway	GPCR downstream signaling; Olfactory Signaling Pathway; Olfactory transduction; Signaling by GPCR
Function	G-protein coupled receptor activity; olfactory receptor activity; receptor activity;