



Anti-OR4P4 (C-terminal) polyclonal antibody (DPABH-01156)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Odorant receptor.
Immunogen	Synthetic peptide derived from C the terminal region of Human OR4P4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituents: 49% PBS, 50% Glycerol, 0.88% Sodium chloride
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [OR4P4 olfactory receptor, family 4, subfamily P, member 4 \[Homo sapiens \]](#)

Official Symbol	OR4P4
Synonyms	OR4P4; olfactory receptor, family 4, subfamily P, member 4; OR4P3P; olfactory receptor 4P4; olfactory receptor 4P3; olfactory receptor, family 4, subfamily P, member 3 pseudogene;
Entrez Gene ID	81300
Protein Refseq	NP_001004124
UniProt ID	Q8NGL7
Chromosome Location	11q11
Pathway	GPCR downstream signaling; Olfactory Signaling Pathway; Olfactory transduction; Signal Transduction; Signaling by GPCR;
Function	G-protein coupled receptor activity; olfactory receptor activity; receptor activity; signal transducer activity;