



Anti-ZNF225 (C-terminal) polyclonal antibody (DPABH-01188)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	Synthetic peptide derived from the C terminal region of Human ZNF225.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituents: 49% PBS, 50% Glycerol, 0.88% Sodium chloride. Note: PBS is without Mg2+ and Ca2+
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C.

GENE INFORMATION

Gene Name [ZNF225 zinc finger protein 225 \[Homo sapiens \]](#)

Official Symbol	ZNF225
Synonyms	ZNF225; zinc finger protein 225; C2-H2 type zinc finger protein; MGC119735;
Entrez Gene ID	7768
Protein Refseq	NP_037494
UniProt ID	Q9UK10
Chromosome Location	19q13.31
Pathway	Gene Expression; Generic Transcription Pathway.
Function	DNA binding; metal ion binding; zinc ion binding;