



Anti-THAP4 (internal region) polyclonal antibody (DPABH-01204)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	THAP4 (THAP domain containing 4) is a protein-coding gene. Diseases associated with THAP4 include tuberculosis. GO annotations related to this gene include DNA binding and metal ion binding.
Immunogen	Synthetic peptide derived from an internal region of Human THAP4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituents: 49% PBS, 50% Glycerol, 0.88% Sodium chloride. Note: PBS is without Mg2+ and Ca2+
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C.

GENE INFORMATION

Gene Name	THAP4 THAP domain containing 4 [Homo sapiens]
Official Symbol	THAP4
Synonyms	THAP4; THAP domain containing 4; THAP domain-containing protein 4; CGI 36; PP238; CGI-36;
Entrez Gene ID	51078
Protein Refseq	NP_001157828
UniProt ID	Q8WY91
Chromosome Location	2q37.3
Function	DNA binding; metal ion binding;