



Anti-TXNRD3NB (aa 95-121) polyclonal antibody (DPABH-01755)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TXNRD3NB (thioredoxin reductase 3 neighbor) is a protein-coding gene.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 95-121 of Human TXNRD3IT1 (NP_001034872.1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	Constituent: 99% PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C or at -20°C for long term storage.

GENE INFORMATION

Gene Name [TXNRD3NB thioredoxin reductase 4 neighbor \[Homo sapiens \]](#)

Official Symbol	TXNRD3NB
Synonyms	TXNRD3NB; thioredoxin reductase 3 neighbor; TR2IT1; TXNRD3IT1; TXNRD3NT1; protein TXNRD3NB; TXNRD3 neighbor gene protein; thioredoxin reductase 3 new transcript 1; thioredoxin reductase 2 intronic transcript 1; thioredoxin reductase 3 intronic transcript 1; thioredoxin reductase 3 neighbor gene protein;
Entrez Gene ID	645840
Protein Refseq	NP_001034872.1
UniProt ID	Q6F5E7
Pathway	Selenium Metabolism and Selenoproteins.