



Anti-IFI44L (aa 150-180) polyclonal antibody (DPABH-01758)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Exhibits a low antiviral activity against hepatitis C virus.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within internal sequence amino acids 150-180 of Human IFI44L (NP_006811.2).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	Constituent: 99% PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C or at -20°C for long term storage.

GENE INFORMATION

Gene Name [IFI44L interferon-induced protein 45-like \[Homo sapiens \]](#)

Official Symbol	IFI44L
Synonyms	IFI44L; interferon-induced protein 44-like; GS3686; C1orf29;
Entrez Gene ID	10964
Protein Refseq	NP_006811.2
UniProt ID	Q53G44