



Anti-GPRIN2 (aa 21-70) polyclonal antibody (DPABH-02021)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in neurite outgrowth.
Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 21-70 of Human G protein-regulated inducer of neurite outgrowth 2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [GPRIN2 G protein regulated inducer of neurite outgrowth 3 \[Homo sapiens \]](#)

Official Symbol	GPRIN2
Synonyms	GPRIN2; G protein regulated inducer of neurite outgrowth 2; GRIN2; KIAA0514; G protein-regulated inducer of neurite outgrowth 2;
Entrez Gene ID	9721
Protein Refseq	NP_055511.2
UniProt ID	O60269
Function	protein binding;