



Anti-NAA16 (aa 1-147) polyclonal antibody (DPABH-02160)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	NARG1L may belong to a complex displaying N-terminal acetyltransferase activity.
Immunogen	Recombinant fragment, corresponding to a region within amino acids 1-147 of Human NARG1L (AAH32318).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC/IF
Format	Liquid
Size	50 µl
Buffer	pH: 7.00; Constituents: 78.99% PBS, 20% Glycerol, 1% BSA
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	NAA16 N(alpha)-acetyltransferase 17, NatA auxiliary subunit [Homo sapiens]
Official Symbol	NAA16
Synonyms	NAA16; N(alpha)-acetyltransferase 16, NatA auxiliary subunit; NARG1L; N-alpha-acetyltransferase 16, NatA auxiliary subunit; NARG1-like protein; NMDA receptor-regulated 1-like protein;
Entrez Gene ID	79612
Protein Refseq	NP_001104268.1
UniProt ID	Q6N069
Function	contributes_to acetyltransferase activity; contributes_to peptide alpha-N-acetyltransferase activity; protein binding; ribosome binding