



Anti-CCDC174 (aa 19-118) polyclonal antibody (DPABH-03055)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment corresponding to amino acids 19-118 of Human C3orf19.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	CCDC174 coiled-coil domain containing 175 [Homo sapiens]
Official Symbol	CCDC174

Synonyms	CCDC174; coiled-coil domain containing 174; C3orf19; coiled-coil domain-containing protein 174;
Entrez Gene ID	51244
Protein Refseq	NP_057558.3
UniProt ID	Q6PII3