



Anti-PET117 (aa 22-74) polyclonal antibody (DPABH-03160)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment, corresponding to amino acids 22-74 of Human PET117 (UniProt Q6UWS5).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	PET117 PET118 homolog (S. cerevisiae) [Homo sapiens]
Official Symbol	PET117

Synonyms	PET117; PET117 homolog (<i>S. cerevisiae</i>); CSRP2BP; protein PET117 homolog, mitochondrial; alternative protein CSRP2BP;
Entrez Gene ID	100303755
Protein Refseq	NP_001158283.1
UniProt ID	L0R6F6