



Anti-RNASE11 (aa 14-199) polyclonal antibody (DPABH-03473)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RNASE11 belongs to the pancreatic ribonuclease family.
Immunogen	Recombinant fragment, corresponding to a region within amino acids 14-199 of Human RNASE11 (UniProt ID: Q8TAA1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.00; Constituents: 78% PBS, 20% Glycerol, 1% BSA
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	RNASE11 ribonuclease, RNase A family, 11 (non-active) [Homo sapiens]
------------------	--

Official Symbol	RNASE11
Synonyms	RNASE11; ribonuclease, RNase A family, 11 (non-active); C14orf6, chromosome 14 open reading frame 6; probable ribonuclease 11; RNase 11; ribonuclease 11; C14orf6;
Entrez Gene ID	122651
Protein Refseq	NP_660293
UniProt ID	Q5GAN5
Chromosome Location	14q11.1
Function	endonuclease activity; hydrolase activity; nucleic acid binding; pancreatic ribonuclease activity;