



Anti-CHD2 (aa 260-476) polyclonal antibody (DPABH-03736)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Sequence-selective DNA-binding protein.
Immunogen	Recombinant fragment corresponding to a region within N terminal amino acids 260-476 of Human CHD2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.00; Constituents: 0.75% Glycine, 1.21% Tris, 10% Glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [CHD2 chromodomain helicase DNA binding protein 2 \[Homo sapiens \]](#)

Official Symbol	CHD2
Synonyms	CHD2; chromodomain helicase DNA binding protein 2; chromodomain-helicase-DNA-binding protein 2; DKFZp547I1315; DKFZp686E01200; DKFZp781D1727; FLJ38614; CHD-2; ATP-dependent helicase CHD2;
Entrez Gene ID	1106
Protein Refseq	NP_001036037
UniProt ID	O14647
Chromosome Location	15q26
Function	ATP binding; ATP-dependent DNA helicase activity; DNA binding; helicase activity; hydrolase activity; nucleotide binding;