



Anti-MLLT10 (aa 281-293) polyclonal antibody (DPABH-03996)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Probably involved in transcriptional regulation. In vitro or as fusion protein with MLL has transactivation activity. Binds to cruciform DNA.
Immunogen	Synthetic peptide (C-RLEDTTARFTNAN) corresponding to internal sequence amino acids 281-293 of Human AF10 (NP_004632.1).
Isotype	IgG
Source/Host	Goat
Species Reactivity	Rat
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	200 µl
Buffer	Constituents: 0.5% BSA, Tris saline, pH 7.3
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	MLLT10 myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila); translocated to, 11 [Homo sapiens]
Official Symbol	MLLT10
Synonyms	MLLT10; myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila); translocated to, 10; AF10; protein AF-10; type I AF10 protein; type IV AF10 protein; type III AF10 protein; ALL1-fused gene from chromosome 10 protein;
Entrez Gene ID	8028
Protein Refseq	NP_001182555.1
UniProt ID	P55197
Function	DNA binding; protein binding; sequence-specific DNA binding transcription factor activity; zinc ion binding