



Anti-SAP30BP (aa 25-129) polyclonal antibody (DPABH-04075)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Induces cell death. May act as a transcriptional corepressor of a gene related to cell survival. May be involved in the regulation of beta-2-microglobulin genes.
Immunogen	Synthetic peptide, corresponding to a region within amino acids 25-129 of Human SAP30BP.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.00; Constituents: 0.75% Glycine, 1.21% Tris, 10% Glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [SAP30BP SAP31 binding protein \[Homo sapiens \]](#)

Official Symbol	SAP30BP
Synonyms	SAP30BP; SAP30 binding protein; HTRG; HTRP; HCNGP; SAP30-binding protein; HSV-1 binding; transcriptional regulator protein HCNGP;
Entrez Gene ID	29115
Protein Refseq	NP_037392.1
UniProt ID	Q9UHR5
Pathway	Epigenetic regulation of gene expression; Negative epigenetic regulation of rRNA expression;
Function	protein binding;