



Anti-KRT6B (aa 134-564) polyclonal antibody (DPABH-04302)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	KRT6B (keratin 6B) is a protein-coding gene. Diseases associated with KRT6B include pachyonychia congenita, and pachyonychia congenita type 2. GO annotations related to this gene include structural constituent of cytoskeleton. An important paralog of this gene is KRT74.
Immunogen	Recombinant fragment corresponding to a region within amino acids 134-564 of Human KRT6B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.00; Constituents: 79% PBS, 20% Glycerol
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	KRT6B keratin 6B [Homo sapiens]
Official Symbol	KRT6B
Synonyms	KRT6B; keratin 6B; keratin like 1 (a type II keratin sequence); KRTL1; keratin, type II cytoskeletal 6B; keratin-6B; cytokeratin 6B; cytokeratin-6B; type-II keratin Kb10; keratin, epidermal, type II, K6B; keratin-like 1 (a type II keratin sequence); K6B; PC2; CK6B; CK-6B; KRTL1;
Entrez Gene ID	3854
Protein Refseq	NP_005546
UniProt ID	P04259
Chromosome Location	12q13.13
Function	structural constituent of cytoskeleton;