



Anti-TMEM66 (aa 1-211) polyclonal antibody (DPABH-04480)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM66 (transmembrane protein 66) is a protein-coding gene. GO annotations related to this gene include protein binding and store-operated calcium channel activity.
Immunogen	Recombinant fragment, corresponding to a region within amino acids 1-211 of Human TMEM66 (Uniprot ID: Q96BY9).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.00; Constituents: 78% PBS, 20% Glycerol, 1% BSA
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	TMEM66 transmembrane protein 66 [Homo sapiens]
Official Symbol	TMEM66
Synonyms	TMEM66; transmembrane protein 66; MGC8721; protein FOAP-7; SARAF long isoform; SARAF short isoform; HBV XAg-transactivated protein 3; HBV X-transactivated gene 3 protein; XTP3; FOAP-7; HSPC035; FLJ22274;
Entrez Gene ID	51669
Protein Refseq	NP_057211
UniProt ID	Q96BY9
Chromosome Location	8p12