



Rabbit Anti-Human MPO Polyclonal Antibody (DPABH-05549)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-310 of human MPO (NP_000241.1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Positive Control	Jurkat,HepG2,HL-60,Mouse pancreas,Rat pancreas
Format	Liquid
Size	50 µl, 100 µl
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction Myeloperoxidase (MPO) is a heme protein synthesized during myeloid differentiation that

constitutes the major component of neutrophil azurophilic granules. Produced as a single chain precursor, myeloperoxidase is subsequently cleaved into a light and heavy chain. The mature myeloperoxidase is a tetramer composed of 2 light chains and 2 heavy chains. This enzyme produces hypohalous acids central to the microbicidal activity of neutrophils.

Keywords	MPO; myeloperoxidase;
-----------------	-----------------------

GENE INFORMATION

Entrez Gene ID	4353
-----------------------	----------------------

UniProt ID	P05164
-------------------	------------------------
