



Anti-ASTN2 (full length) polyclonal antibody (DPABH-06509)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ASTN2 (astrotactin 2) is a protein-coding gene. Diseases associated with ASTN2 include migraine without aura, and schizophrenia. An important paralog of this gene is ASTN1.
Immunogen	Full length protein, corresponding to amino acids 1-440 of Human ASTN2 (AAH29272.1).
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	pH: 7.20; Constituent: 99% PBS
Preservative	None
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	ASTN2 astrotactin 2 [Homo sapiens]
Official Symbol	ASTN2
Synonyms	ASTN2; astrotactin 2; astrotactin-2; KIAA0634; bA264C15.1; bA67K19.1;
Entrez Gene ID	23245
Protein Refseq	NP_001171663
UniProt ID	B7ZKP4
Chromosome Location	9q33