



Anti-ZNF436 (full length) polyclonal antibody (DPABH-07671)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be a transcriptional repressor.
Immunogen	Full length protein, corresponding to amino acids 1-470 of Human ZNF436 (NP_085137.1).
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, ICC/IF
Format	Liquid
Size	50 µg
Buffer	pH: 7.20; Constituent: 99% PBS
Preservative	None
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [ZNF436 zinc finger protein 436 \[Homo sapiens \]](#)

Official Symbol	ZNF436
Synonyms	ZNF436; zinc finger protein 436; KIAA1710; Zfp46; DNA-binding protein; ZNF;
Entrez Gene ID	80818
Protein Refseq	NP_001070663
UniProt ID	A0A024RAF3
Chromosome Location	1p36
Pathway	Gene Expression; Generic Transcription Pathway;
Function	DNA binding; metal ion binding; nucleic acid binding; zinc ion binding;