



Anti-GPRIN2 (full length) polyclonal antibody (DPABH-07971)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in neurite outgrowth.
Immunogen	Full length protein, corresponding to amino acids 1-458 of Human G protein-regulated inducer of neurite outgrowth 2 (AAH11672.1).
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µl
Preservative	None
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [GPRIN2 G protein regulated inducer of neurite outgrowth 3 \[Homo sapiens \]](#)

Official Symbol	GPRIN2
Synonyms	GPRIN2; G protein regulated inducer of neurite outgrowth 2; GRIN2; KIAA0514; G protein-regulated inducer of neurite outgrowth 2;
Entrez Gene ID	9721
Protein Refseq	NP_055511.2
UniProt ID	O60269
Function	protein binding;