



## Anti-SELV (aa 303-332) polyclonal antibody (DPABH-08539)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide within Human SELV aa 303-332 (C terminal) conjugated to Keyhole Limpet Haemocyanin (KLH). The exact sequence is proprietary. Database link: P59797
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ICC/IF, Flow Cyt, IHC-P
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	Constituent: 99% PBS
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Store at 4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

### BACKGROUND

**Introduction** May be involved in a redox-related process.

**Keywords**

SELV; selenoprotein V; selenoprotein V;

---

## GENE INFORMATION

**Entrez Gene ID**

[348303](#)

---

**UniProt ID**

[P59797](#)

---