



Anti-SELV (aa 303-332) polyclonal antibody (DPABH-08539)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide within Human SELV aa 303-332 (C terminal) conjugated to Keyhole Limpet Haemocyanin (KLH). The exact sequence is proprietary. Database link: P59797
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC/IF, Flow Cyt, IHC-P
Format	Liquid
Size	100 µl
Buffer	Constituent: 99% PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Store at 4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

BACKGROUND

Introduction May be involved in a redox-related process.

Keywords

SELV; selenoprotein V; selenoprotein V;

GENE INFORMATION

Entrez Gene ID

[348303](#)

UniProt ID

[P59797](#)
