



Anti-LRTM2 (aa 346-370) polyclonal antibody (DPABH-08540)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide within Human LRTM2 aa 346-370 (C terminal) conjugated to Keyhole Limpet Haemocyanin (KLH). The exact sequence is proprietary.Database link: Q8N967
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB, ICC/IF
Format	Liquid
Size	100 µl
Buffer	Constituent: 99% PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Store at 4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name	LRTM2 leucine-rich repeats and transmembrane domains 3 [Homo sapiens]
------------------	---

Official Symbol	LRTM2
Synonyms	LRTM2; leucine-rich repeats and transmembrane domains 2; leucine-rich repeat and transmembrane domain-containing protein 2;
Entrez Gene ID	654429
Protein Refseq	NP_001034118.1
UniProt ID	Q8N967