



Anti-ZCCHC8 (aa 559-590) polyclonal antibody (DPABH-10372)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in pre-mRNA splicing.
Immunogen	Synthetic peptide within Human ZCCHC8 aa 559-590 (C terminal) conjugated to Keyhole Limpet Haemocyanin (KLH). The exact sequence is proprietary. Database link: Q6NZY4
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	Constituent: 99% PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Store at 4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name	ZCCHC8 zinc finger, CCHC domain containing 8 [Homo sapiens]
Official Symbol	ZCCHC8
Synonyms	ZCCHC8; zinc finger, CCHC domain containing 8; zinc finger CCHC domain-containing protein 8; DKFZp434E2220;
Entrez Gene ID	55596
Protein Refseq	NP_060082
UniProt ID	Q6NZY4
Chromosome Location	12q24.31
Function	metal ion binding; nucleic acid binding; protein binding; zinc ion binding;