



Anti-PLEKHA6 (aa 709-736) polyclonal antibody (DPABH-10498)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide within Human PEPP3 aa 709-736 (internal sequence) conjugated to Keyhole Limpet Haemocyanin (KLH). The exact sequence is proprietary. Database link: Q9Y2H5
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	Constituent: 99% PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Store at 4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name	PLEKHA6 pleckstrin homology domain containing, family A member 7 [Homo sapiens]
------------------	---------------------------------------------------------------------------------------------------

Official Symbol	PLEKHA6
Synonyms	PLEKHA6; pleckstrin homology domain containing, family A member 6; PEPP3; PEPP-3; pleckstrin homology domain-containing family A member 6; PH domain-containing family A member 6; phosphoinositol 3-phosphate-binding protein 3; phosphoinositol 3-phosphate-binding protein-3;
Entrez Gene ID	22874
Protein Refseq	NP_055750.2
UniProt ID	Q9Y2H5