



Anti-TRIM74 (aa 165-193) polyclonal antibody (DPABH-10621)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide within Human TRIM50C aa 165-193 (C terminal) conjugated to Keyhole Limpet Haemocyanin (KLH). The exact sequence is proprietary. Database link: Q86UV6
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	Constituent: 99% PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Store at 4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name [TRIM74 tripartite motif containing 75 \[Homo sapiens \]](#)

Official Symbol	TRIM74
Synonyms	TRIM74; tripartite motif containing 74; TRIM50C; tripartite motif-containing protein 74; tripartite motif-containing 74; tripartite motif-containing 50C; tripartite motif-containing protein 50C;
Entrez Gene ID	378108
Protein Refseq	NP_942150.1
UniProt ID	Q86UV6
Function	zinc ion binding;