



Anti-TBC1D2B (aa 570-619) polyclonal antibody (DPABH-12714)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May act as a GTPase-activating protein.
Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 570-619 (EDALQVESQE QPEQAFVKPH LVSEYDIYGF RTVPEDDEEE KLVAKVRALD) of Human TBC1D2B (NP_055894).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 97% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	TBC1D2B TBC1 domain family, member 3B [Homo sapiens]
Official Symbol	TBC1D2B
Synonyms	TBC1D2B; TBC1 domain family, member 2B; TBC1 domain family member 2B;
Entrez Gene ID	23102
Protein Refseq	NP_055894.6
UniProt ID	Q9UPU7
Function	Rab GTPase activator activity; protein binding;