



Anti-SCUBE1 (aa 63-122) polyclonal antibody (DPABH-12769)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Could function as an adhesive molecule and its matrix bound and soluble fragments may play a critical role in vascular biology.
Immunogen	Synthetic peptide within Human SCUBE1 aa 63-122 (N terminal). The exact sequence is proprietary. (NP_766638)Sequence: PGYKGEKQCEDIDECENDYYNGGCVHECINIPGNYRCTCFDGFMLAHDG Database link: Q8IWY4
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Preservative	None
Storage	Shipped at 4°C. Store at 4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name	SCUBE1 signal peptide, CUB domain, EGF-like 1 [Homo sapiens]
Official Symbol	SCUBE1
Synonyms	SCUBE1; signal peptide, CUB domain, EGF-like 1; signal peptide, CUB and EGF-like domain-containing protein 1; signal peptide-CUB domain-EGF-related 1;
Entrez Gene ID	80274
Protein Refseq	NP_766638
UniProt ID	A0JP65
Chromosome Location	22q13
Function	calcium ion binding; identical protein binding; protein heterodimerization activity;