



Rabbit Anti-Human POSTN Polyclonal Antibody (DPABH-12789)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Polyclonal to POSTN
Target	POSTN
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	This antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Conjugate	Unconjugated
Applications	WB, ELISA
Molecular Weight	91 kDa
Format	Liquid
Concentration	Lot specific
Size	100 µl, 200 µl
Buffer	PBS containing 50% glycerol
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, and avoid repeat freeze-thaw cycles.

BACKGROUND

Introduction	Binds to heparin. Induces cell attachment and spreading and plays a role in cell adhesion. May play a role in extracellular matrix mineralization.
Keywords	POSTN;periostin, osteoblast specific factor;PN;OSF2;OSF-2;PDLPOSTN;periostin;periostin;periodontal ligament-specific periostin;osteoblast specific factor 2 (fasciclin I-like);

GENE INFORMATION

Entrez Gene ID	10631
UniProt ID	Q15063