



Anti-FBXL4 (aa 552-565) polyclonal antibody (DPABH-13769)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	FBXL4 (F-box and leucine-rich repeat protein 4) is a protein-coding gene. An important paralog of this gene is FBXL2.
Specificity	DPABH-13769 shows no cross reactivity with other proteins.
Immunogen	Synthetic peptide corresponding to Human FbxL4 aa 552-565 (C terminal).Sequence: DELASNCTRLQQLD Database link: Q9UKA2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µg
Buffer	Constituents: 2.5% BSA, 0.45% Sodium chloride, 0.1% Dibasic monohydrogen sodium phosphate
Preservative	0.025% Sodium Azide
Storage	Shipped at 4°C. Store at 4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name	FBXL4 F-box and leucine-rich repeat protein 4 [Homo sapiens]
Official Symbol	FbxL4
Synonyms	FBXL4; F-box and leucine-rich repeat protein 4; F-box/LRR-repeat protein 4; FBL4; FBL5; F-box protein FBL4; F-box protein FBL4/FBL5;
Entrez Gene ID	26235
Protein Refseq	NP_036292
UniProt ID	Q9UKA2
Chromosome Location	6q16.1-q16.3