



Anti-ZNF407 (aa 2000-2050) polyclonal antibody (DPABH-15048)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	Synthetic peptide within Human ZNF407 aa 2000-2050. The exact sequence is proprietary. NP_060227.2Sequence: GEVEGRAGLEEQGRPGAKDVLIQLPGQEVSHVAADPEAPEIQMFPQAQES P Database link: Q9C0G0
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IP
Format	Liquid
Size	100 µg
Buffer	Constituent: 99% Tris citrate/phosphate, pH 7 - 8
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Store at 4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.

GENE INFORMATION

Gene Name	ZNF407 zinc finger protein 407 [Homo sapiens]
Official Symbol	ZNF407
Synonyms	ZNF407; zinc finger protein 407; FLJ13839; FLJ20307; KIAA1703;
Entrez Gene ID	55628
Protein Refseq	NP_001139661
UniProt ID	Q9C0G0
Chromosome Location	18q23
Function	DNA binding; metal ion binding; zinc ion binding;