



Anti-PRDM14 (aa 217-266) polyclonal antibody (DPABH-15441)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	A region within a synthetic peptide: ASLHHAISGL LVPPDSSGSD SLPQTLDKDS LQLPEGLCLM QTVFGEVPHF, corresponding to internal sequence amino acids 217-266 of Human PRDM14
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC-P
Format	Liquid
Size	100 µg
Buffer	Constituents: 2% Sucrose, PBS.
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	PRDM14 PR domain containing 14 [Homo sapiens]
Official Symbol	PRDM14
Synonyms	PRDM14; PR domain containing 14; PR domain zinc finger protein 14; PR domain-containing protein 14; PFM11; MGC59730;
Entrez Gene ID	63978
Protein Refseq	NP_078780
UniProt ID	Q9GZV8
Chromosome Location	8q13.3
Function	DNA binding; metal ion binding; zinc ion binding;